



PHILIPS

Customer Service
Solutions



Simplify your ultrasound quality control process

The Philips Ultrasound Accreditation Assist Program



In today's challenging healthcare environment, maintaining the quality and performance of your medical imaging systems is critical. The American College of Radiology (ACR) and other accrediting bodies require healthcare facilities to create a standardized quality assurance program for imaging and patient safety. The Philips team of highly skilled Clinical Education Specialists and Field Service Engineers can help you meet your accreditation needs with our comprehensive Ultrasound Accreditation Assist Program. We'll help simplify your ultrasound quality control process and maintain your accreditation status. The benefit for you is more time to focus on patient care.

Program Overview

- Skilled Philips Field Service Engineers will perform an annual system performance evaluation meeting all ACR requirements in conjunction with a scheduled Planned Maintenance service. Eliminating a separate meeting allows you more time for patient care.*
- Professionally trained Philips Clinical Education Specialists will complete the required quality control and provide consultation on accreditation processes and documentation. Based on your staffing capacity, they can also train your sonographers to complete the testing.
- As a bonus, with the Ultrasound Accreditation Assist Program, you will receive login access to the introductory e-learning course “Quality Control in Ultrasound Lab Accreditation”.
- With a Clinical Education Flex Account, you can reserve funds for continuing education for your staff accreditation needs.



Program Options

	Service <i>We take care of Routine QC</i> 	Training <i>We prepare your staff to perform Routine QC</i> 
Routine Quality Control	Performed annually by Philips Clinical Education Specialists	Your staff will learn how to perform testing in a 1-day onsite training session
Planned Maintenance / Annual Survey	Entitled in service agreement (if any) performed by Philips Field Service Engineers	Entitled in service agreement (if any) performed by Philips Field Service Engineers
Accreditation Consultation	Included	Included
E-learning Course	Included	Included
Clinical Education Flex Account	Optional	Optional

*Planned Maintenance is an entitlement of a service agreement. If no service agreement is in place, additional Planned Maintenance is required.

For more information: Please call **1 800 522 7022**, or visit the Philips Ultrasound Accreditation Assist web page at www.learningconnection.philips.com/ultrasound-accreditation-assist

© 2016 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

www.philips.com



Please visit www.usa.philips.com/healthcare-medical-education

APR 2016