



# 3D Auto LAA eLearning

3D Auto LAA eLearning – Philips 1st dedicated solution for live LAA sizing. The measurements provided are Perimeter, Area and Max and min diameters.

This eLearning will review the features and workflow of the 3D Auto LAA tool.

The 3D Auto LAA eLearning module will:

- Explain what Atrial Fibrillation is and its relationship to Stroke.
- Describe the 3D Auto Left Atrial Appendage (LAA) tool
- What are the clinical benefits of using this tool
- Demonstrate how to measure the LAA with the 3D Auto LAA tool
- Describe the editing options within the tool

This course is intended for customers who have purchased the 3D Auto LA Aoption.

This ASRT-accredited module is a self-paced tutorial consisting of enlightening content and colorful images.



### For more information

Contact a Philips ultrasound clinical services coordinator at 800-522-7022 and visit our education catalog at

[www.learningconnection.philips.com/ultrasound](http://www.learningconnection.philips.com/ultrasound)